



BC Cancer Agency

CARE & RESEARCH

An agency of the Provincial Health Services Authority



CANADA'S MICHAEL SMITH
**GENOME
SCIENCES**
CENTRE

Detecting structural variants

February 14th 2018

Karen Mungall, Bioinformatics Coordinator - Analysis Groups

Canada's Michael Smith Genome Sciences Centre

BC Cancer Agency

- Part 1
 - What is a gene fusion
 - How do they arise
 - Why are they important
- Part 2
 - Considerations for tool selection
 - What to do with the data
- Part3
 - Comprehensive SV detection



BC Cancer Agency

CARE & RESEARCH


An agency of the Provincial Health Services Authority

PART 1


PART 1

Gene fusions - What, how, why?


- When two separate genes come together to form a new chimeric gene. The resulting protein product may lead to abnormal expression levels and function and may in turn cause the abnormal proliferation of cells and cancer development



How does a gene fusion arise?




The first fusion gene identified, is known as the Philadelphia chromosome. It arises from a translocation event involving the 5' part of the BCR gene on chr22 fusing to the 3' part of the ABL1 gene on chr9. It was first discovered in chronic myelogenous leukemia (CML). *BCR-ABL1* has been found to occur in more than 95% of CML patients and to exert its oncogenic phenotype by encoding a constitutively active ABL1 kinase






Gene fusion time line






Recurrency in gene fusions



Trends in fusion functionality


A Gene fusion landscapes are diverse

The diversity, abundance, and connection to etiology of gene fusions varies across both cancers and individuals




B Gene fusion networks elucidate fusion pairings

Network studies show that most fusion genes fuse with very few partners, and that different cancer types have signature fusion networks




C The frequency of fusions in cancers varies considerably

Fusions tend to be rare, but can be predominant, and anti-correlate with other somatic mutations



D Fusion genes tend to have specific functions


Molecular functions relating to kinase or DNA-binding activity are enriched in genes forming fusions



Structural features of fusion proteins


A Breakpoint locations tend to preserve protein function

Breakpoints tend to occur in disordered regions and maintain reading frames and protein globularity




B Fusion proteins are relatively depleted in domains

Proteins which form fusions have fewer domains than other proteins, but fusion transcripts encode more domains than expected by chance.




C Fusion proteins contain specific domain architectures

Domain recombinations in fusion proteins are non-random and sometimes novel



D Disorder may contribute to fusion protein functionality

The increased disorder in fusion proteins could promote the viable joining of different domains and offer flexibility for internal interactions



- **TMPRSS2-ERG Fusion Gene in Prostate Cancer**

High expression of TMPRSS2-ERG gene fusion together with prostate-specific antigen levels are indicators for likelihood of recurrence and shortened time to recurrence

Prognostic Significance of *TMPRSS2-ERG* Fusion Gene in Prostate Cancer ⇒





Figure 2.

Relation of TMPRSS2-ERG fusion transcript expression in combination with preoperative serum PSA level to time to recurrence (Kaplan-Meier curves). A combination of high PSA level and high TMPRSS2-ERG expression was associated with the shortest time to recurrence.

- Confirmation of diagnosis
 - BCR-ABL1 in CML patients hall mark fusion seen in ~95% patients
- Specific subgroup:
 - EML4-ALK fusion is seen in around 5% of NSCLC patients





- Patients with *ALK* rearrangements **do not** benefit from EGFR-specific TKI therapy but may be considered for therapy targeting the constitutively activated receptor tyrosine kinase that results from *EML4-ALK* and other *ALK* fusions. Crizotinib is the first FDA-approved ALK TKI. It is indicated for treatment of locally advanced or metastatic NSCLC in patients whose tumors are positive for ALK as determined using an FDA-approved test.
- Additionally, *EGFR*, *KRAS*, and *ALK* mutations are almost always mutually exclusive (ie, mutations of only 1 of the 3 genes occur within any individual tumor).
- Methods for detecting the *ALK* rearrangements include FISH, PCR, and immunohistochemical (IHC) staining. ALK tests are often run in conjunction with tests for EGFR and KRAS mutations
- Outcome: Sensitive to ALK inhibitors eg Crizotinib Resistant to EGFR Tyrosine Kinase Inhibitors



BC Cancer Agency

CARE & RESEARCH

An agency of the Provincial Health Services Authority


PART 2

PART 2

Tool selection and what to do with results


How to detect fusions

Fluorescent in situ Hybridization (FISH)




How to detect fusions

Polymerase Chain Reaction (PCR)



How to detect fusions

Massively parallel sequencing (Illumina sequencing by synthesis)




- Advantages
 - specific and fast
- Disadvantages
 - need to know ahead of time what you want to find



Identify actionable event
(can be SNV, indel, fusion breakpoint,
fusion gene partners or exon in gene)




Create probe sequence containing event



Screen probe sequence with target fastq
or bam file


– Assembly based fusion detection

- Advantages
 - Comprehensive event detection
 - Higher specificity
 - generates contigs for better interrogation of event breakpoint
- Disadvantages
 - Large resource requirement with multiple steps and slow
 - Lower sensitivity



– Alignment based fusion detection



- Advantages
 - Fast with lower resource requirement
 - Higher sensitivity
- Disadvantages
 - lower specificity






- Sample source
 - Fresh or FFPE
- Data type
 - DNA or RNA
- Input data
 - Fastq or bam
 - readlength
- Test samples
 - Individual matched samples eg DNA, RNA
 - Multiple individual analysis eg trio

-
- Sensitivity
 - Specificity
 - Speed
 - Resources
 - Deterministic




Genome



Transcriptome

Gene fusion visualization





-
- <http://www.tumorfusions.org/>
 - COSMIC
 - Mitleman



TUMOR FUSION GENE DATA PORTAL

Landscape of cancer-associated fusions using the Pipeline for RNA sequencing Data Analysis

Transcripts fusion as a result of genomic rearrangement is an important class of somatic alteration, as a cancer initiating event and as a molecular therapeutic target for specific tumors. Our [Pipeline for RNA sequencing Data Analysis \(PRADA\)](#) enables us to detect fusion transcripts with high confidence comprehensively. Based on integrated analysis of paired-end RNA sequencing and DNA copy number data from [The Cancer Genome Atlas \(TCGA\)](#), The Tumor Fusion Gene Data Portal provides a bona-fide fusion list across many tumor types.




Figure 1. Fusion transcripts. Fusion transcripts are chimeric mRNAs encoded from the joined parts of two genes, and may occur as a result of genomic rearrangements.





Figure 2. Detection of fusion transcripts. PRADA detects fusion transcripts through identification of discordant read pairs and junction spanning reads.

Navigation
Introduction
Summary
By Gene
By Fusion
By Sample
Cancer Type
Normal Tissue
Download
Help
Contact Us



The screenshot shows the COSMIC website homepage. At the top is the COSMIC logo and navigation menu. Below the menu is a search bar with a search button. The main content area features a heading for the latest release, a brief description of COSMIC, and a search example. Below this is a 'Projects' section with four items: COSMIC, Cell Lines Project, COSMIC-3D, and Cancer Gene Census. At the bottom is a 'Data curation' section with four links.

COSMIC
Catalogue Of Somatic Mutations In Cancer

Projects ▾ Data ▾ Tools ▾ News ▾ Help ▾ About ▾ Genome Version ▾ Search COSMIC... **SEARCH**

COSMIC v83, released 07-NOV-17

COSMIC, the Catalogue Of Somatic Mutations In Cancer, is the world's largest and most comprehensive resource for exploring the impact of somatic mutations in human cancer.

Start using COSMIC by searching for a gene, cancer type, mutation, etc. below.

eg Braf, COLO-829, Carcinoma, V600E, BRCA-UK, Campbell **SEARCH**

Projects

COSMIC is divided into several distinct projects, each presenting a separate dataset or view of our data:

-  **COSMIC**
The core of COSMIC, an expert-curated database of somatic mutations
-  **Cell Lines Project**
Mutation profiles of over 1,000 cell lines used in cancer research
-  **COSMIC-3D**
An interactive view of cancer mutations in the context of 3D structures
-  **Cancer Gene Census**
A catalogue of genes with mutations that are causally implicated in cancer

Data curation

-  [Gene Curation](#) — details of our manual curation process
-  [Gene Fusion Curation](#) — details of our curation process for gene fusions
-  [Genome Annotation](#) — information on the annotation of genomes
-  [Drug Resistance](#) — curation of mutations conferring drug resistance



NATIONAL CANCER INSTITUTE
U.S. National Institutes of Health | www.cancer.gov

CANCER GENOME ANATOMY PROJECT

CGAP How To
Genes
Chromosomes
Tissues
SAGE Genie
RNAi
Pathways
Tools

Chromosomes

Tools

- [FISH-mapped BACs](#)
- [Genetic and Physical SNP Maps](#)
- [Mitelman Searchers](#)
- [SNP500Cancer](#)

CGAP Data

- [Download](#)

Purchase CGAP Reagents

- [CCAP BAC Clones](#)

Related Links

- [Atlas of Genetics & Cytogenetics in Oncology & Haematology](#)
- [Progenetix](#)
- [SKY/M-FISH & CGH Database](#)

Quick Links:

- [ICG](#)
- [NCI Home](#)
- [NCICB Home](#)
- [NCBI Home](#)
- [OCG](#)

Mitelman Database of Chromosome Aberrations and Gene Fusions in Cancer

Searching the Database

The information in the Mitelman Database of Chromosome Aberrations and Gene Fusions in Cancer relates chromosomal aberrations to tumor characteristics, based either on individual cases or associations. All the data have been manually culled from the literature by Felix Mitelman, Bertil Johansson, and Fredrik Mertens.

CGAP has developed six web search tools to help you analyze the information within the Mitelman Database:

- The **Cases Quick Searcher** allows you to query the individual patient cases using the four major fields: aberration, breakpoint, morphology, and topography.
- The **Cases Full Searcher** permits a more detailed search of the same individual patient cases as above, by including more cytogenetic field choices and adding search fields for patient characteristics and references.
- The **Molecular Biology Associations Searcher** does not search any of the individual patient cases. It searches studies pertaining to gene rearrangements as a consequence of cytogenetic aberrations.
- The **Clinical Associations Searcher** does not search any of the individual patient cases. It searches studies pertaining to clinical associations of cytogenetic aberrations and/or gene rearrangements.
- The **Recurrent Chromosome Aberrations Searcher** provides a way to search for structural and numerical abnormalities that are recurrent, i.e., present in two or more cases with the same morphology and topography.
- The **Reference Searcher** queries only the references themselves, i.e., the references from the individual cases and the molecular biology and clinical associations.

Database last updated on December 14, 2017
 Total number of cases = **68,170**
 Total number of gene fusions = **11,124**

Need help! To learn about the Mitelman Database and how to search it, please visit:

- **All about the Mitelman Database**, which provides background information about the Mitelman Database of Chromosome Aberrations and Gene Fusions in Cancer.
- **Mitelman Database Search Help**, which contains information on how to use the search tools.
- **ISCN Abbreviated Terms and Symbols**, which provides a list of terms and symbols used to describe chromosome abnormalities.

Citation of the Database



BC Cancer Agency

CARE & RESEARCH

An agency of the Provincial Health Services Authority

PART 3

PART 3

Comprehensive structural variant detection

Comprehensive structural variant detection

- Multiple tool input
- Clustering of breakpoints
- Consistent breakpoint calling
- Data pairing
- Evidence support
- Standard annotation
- Standard output format




BC Cancer Agency

CARE & RESEARCH

An agency of the Provincial Health Services Authority

MAVIS

Merging, Annotation, Validation, and
Illustration of Structural variants



1. Cluster
2. Filter based on proximity to annotations
 - Call a new merged breakpoint pair from the group
3. Gather read evidence
4. Call breakpoint pairs (contig, split read, flanking pairs)
 - contig, split read, flanking pairs
5. Annotate with gene and transcript level information
 - Build Fusion Transcripts for exact calls
6. Draw SVGs for all calls
7. Pair calls between libraries
 - Somatic, Expressed
8. Summary
 - Standard output file with HGVS nomenclature



Merging

- Takes inputs from any SV caller as long as it is put in a common format
- filters based on user-defined masked regions
- Splits calls by type
- Merges based on proximity
 - uses a clique finding algorithm
 - followed by hierarchical clustering for larger clusters than cannot be computed exactly inexpensively



- Uses bam files to collect support for the input event calls
 - Uses read pair fragment distribution to define intervals of where reads will be collected from
 - collects spanning, split, and half-mapped reads
 - collects flanking and compatible-flanking pairs
 - standardizes cigar/read-alignments to ensure reproducible calls
 - uses a collapsed annotation model to adjust these intervals and calculate read-pair fragment sizes for transcriptomes
- Local assembly
 - Does not attempt to resolve or assemble repeats longer than the read/kmer length
- Calling breakpoint pairs
 - call by contig
 - call by split reads when calling by contig fails
 - call by flanking (or split and flanking) pairs when calling by split reads fails

Gene level

Nearby genes

Genes encompassed by the event


Genes at the breakpoints

Fusion Transcript level

Exon/intron of breakpoint

Uses a splicing model to predict if the fusion will be in or out of frame

Predicts domain retention by re-aligning protein domain sequences to the new amino acid sequence of the fusion protein



For any fusions with breakpoint level resolution a figure and putative splicing products are produced







Illustration- alternative splicing





Want more information?

-
- [MAVIS](#)
 - <https://github.com/bcgsc/mavis/>
 - mavis@bcgsc.ca
 - Submitted to Bioinformatics
 - Poster presentation AGBT 2018



Acknowledgements

Production Bioinformatics Analysis Groups

- Amir Zadeh
- Morgan Bye
- Simon Chan

Production Bioinformatics Database Groups

- Eric Chuah
- Tina Wong
- Dean Cheng
- Kane Tse

Yussanne Ma
Steven Jones
Marco Marra