




BC Cancer Agency

CARE & RESEARCH

An agency of the Provincial Health Services Authority




CANADA'S MICHAEL SMITH
**GENOME
SCIENCES**
CENTRE

Quality Control of Next Generation Sequence Data

January 17, 2018

Kane Tse, Assistant Bioinformatics Coordinator
Canada's Michael Smith Genome Sciences Centre
BC Cancer Agency

Canada's Michael Smith Genome Sciences Centre





CANADA'S MICHAEL SMITH
G E N O M E
SCIENCES
C E N T R E



Genomics & Bioinformatics Research Centre

- Part of the Cancer Research Centre of the BC Cancer Agency





Sequencing Platforms

- Illumina sequence-by-synthesis instruments
 - NextSeq, MiSeq, HiSeq 2500, HiSeqX instruments
- Sanger capillary-based sequencing
 - Life 3730 XL
- Monthly
 - 1,500 libraries constructed
 - >80 terabases sequenced

Bioinformatic Analysis

- 3 large-scale compute clusters
 - 800 nodes, 24,000 hyperthreaded cores, 120TB RAM
- Multiple team-specific clusters
 - Ex - BioQC team: 320 cores, 2.5TB RAM
- 20 Petabytes of storage

Overview





Description

- What is Quality Control?
- How is Quality Control performed?
- Why is Quality Control important to you?


Value


- Quality Control ensures accurate results
- Quality Control can enhance interpretation of results
- Quality Control has scientific merit in publications

Examples

- Unusual cases encountered, and their impact on QC

What is Quality Control





If you don't have time to do it right you must have time to do it over. [Unknown]

BCGSC spends time & effort ensuring Quality


- Many teams monitor quality
 - Tend to be manual checks
 - Relies on experience & expertise for detection
- Bioinformatics Quality Control group
 - Automated pipeline to monitor quality and report issues

Why do we care about Quality?

- Identify potential issues before data analysis begins
- Inform collaborators about their experiment
- Improve our laboratory & bioinformatics processes



Definition



qual-i·ty con·trol

/'kwälədē kən'trōl/

a system of maintaining standards in **manufactured products** by testing a sample of the output against the **specification**.

<http://whatis.techtarget.com/definition/quality-control-QC>

- **manufactured products** = NGS sequence data
- **specification** = type of experiment (WGS, Capture, miRNA)




BC Cancer Agency

CARE & RESEARCH

An agency of the Provincial Health Services Authority

Quality Assurance (QA) vs Quality Control (QC)




CANADA'S MICHAEL SMITH
**GENOME
SCIENCES**
CENTRE



Quality Assurance

- Also plays a big role at the BCGSC
 - But not the focus of today's discussion




Different Levels of Quality Control

- Level 0: Non-Alignment based metrics
- Level 1: Alignment against a reference genome
- Level 2: Assessment after bioinformatic analysis
 - eg. Variant calling, expression quantification




Levels of Quality Control



	Industry Definitions	Bioinformatic Context
Level 0	<ul style="list-style-type: none">• Raw unprocessed data• Directly observed on the instrument• Absolute measurements	<ul style="list-style-type: none">• Input = fastq files from sequencer• Indifferent to protocols, regardless of pipeline (WGS, RNAseq, etc.)
Level 1	<ul style="list-style-type: none">• Quality Controlled data• Associated with metadata• Compared with calibrations	<ul style="list-style-type: none">• Using aligner (BWA or Novoalign) to compare against “standards” (human, mouse reference genomes, etc.)• Mapping rate, dup-rate, paired
Level 2	<ul style="list-style-type: none">• Derived products that require scientific & technical interpretation• Standards defined by the community that collects or utilizes the data	<ul style="list-style-type: none">• Assembly• Expression levels• Variant calling• On-Target Rate

How is Quality Control performed?





Laboratory QC

- DNA Quantification, Agilent traces, Cluster density, intensity, focus scores, PF rate, Q30/Q20, index splitting

Bioinformatic Level 0 QC


- 60 metrics
 - total_reads, contamination, reagent_leftover, miRNA_adapter...

Bioinformatic Level 1 QC

- Alignment (3):
 - % aligned to genome, % properly paired reads, % duplicate rate...
- ChIP-seq (6):
 - Fraction of reads in peaks (FRiP), domain reads as % of mapped reads..
- Bisulfite-seq (4):
 - Lambda bisulfite conversion rate, human bisulfite conversion rate...
- RNAseq (10)
 - Num Genes Covered @ 1X/10X, Percent reads mitochondrial, intergenic reads...
- miRNA (2)
 - Num. miRNA reads, Diversity of miRNA species

QC in the Lab





Pre-Sequencing

- DNA quantification
 - Determine how much DNA is in a sample
- qPCR
 - Determine how many fragments contain Illumina adapters

On Instrument


- First base report
 - Try to detect library issues or machine issues
 - Look for biased libraries from basecalls
 - Review cluster density
- Post-run QC
 - Q30/Q20 scores – contamination of cleavage mix, temperature of instrument
 - Index splitting – uneven pooling, unknown indices

Level 0 QC





Why level 0?



It's Fast

- QC all lanes within 24 hours of sequencing
- Rapid feedback to the lab on go/no-go for subsequent lanes

It's Universal

- Works regardless of protocol or sequencing method
- Detects reagents, spike-ins
- Scan & optionally remove microbial genomes


It's Consistent

- Metrics are generated and loaded automatically into a DB
- Forms a basis for historical comparison & trend analysis




Every lane analyzed for a standard set of metrics

- Some metrics used for pass/fail assessment
- All metrics stored in a database for historical comparison





What can you look for without alignment?



Reagent content

- Detect sequences that contain adapters, vectors, standards, ladders

Microbial Contamination


- Use read classification tools like BioBloomTool (BBT) to detect specific microbial contaminants (45 species)

Index splitting & Pooling problems

- Check if the index no-match bin contains a large number of reads
- Check for expected indexes that are missing reads

Sample Swap

- Compare variant calls between samples of same individual
- Look for spike-ins (PhiX or a GSC-specific spike-in)
- Check that the distribution of indices matches what was pooled



Multiple methods of detection

- SNP Concordance (human libraries only)
- Customized spike ins (WGS, RNAseq, amplicon, WGBS)
- Index splitting (for pooled libraries)

SNP concordance

- Bioinformatic implementation of Affy's 500k chip array

Patient SNP Comparison Table						
	HFJCMCCXY 8 CTAAGG-TATCGCAG	HFJCMCCXY 8 GATATA-AGATCTCG	HFJCMCCXY 8 CTAAGG-TCGACGTA	HFJCMCCXY 8 CTAAGG-ATGATCGA	HFJCMCCXY 8 CTAAGG-GACTTAGC	
	P02636	P02633	P02636	P02636	P02636	
HFJCMCCXY 8 CTAAGG-GACTTAGC	P02636 0.88	P02633 0.657	P02636 0.88	P02636 0.88	P02636 1.0	
HFJCMCCXY 8 CTAAGG-ATGATCGA	P02636 0.885	P02633 0.664	P02636 0.885	P02636 1.0		
HFJCMCCXY 8 CTAAGG-TCGACGTA	P02636 0.888	P02633 0.654	P02636 1.0			
HFJCMCCXY 8 GATATA-AGATCTCG	P02633 0.658	P02633 1.0				
HFJCMCCXY 8 CTAAGG-TATCGCAG	P02636 1.0					

Regenerate Snp Tables

Spike Ins

- Add 200bp oligos into each sample at tiny amounts
- Detect those oligos in sequenced data (~10,000 reads)



Categories of QC metrics



Sequencing Quality

- Adapters, reagents, dimers
- Duplicate rate
- Contamination
- Read quality
- PF Rate (Chastity Passed)
- Coverage

Success of Laboratory Processes

- Bisulfite conversion rate
- Pooling efficiency
- ChIP capture efficiency
- On-target read rate (specific capture)
- Mitochondrial or rRNA content

Sample Degradation

- RNA degradation
- Fragment size

Historical Comparison


- Lane to lane comparison

Sample Identity

- Plasmid Spike Ins
- SNP concordance
- Index splitting

Gene Complexity and Library Diversity

- miRNA diversity
- # of Genes detected
- Intergenic content
- Intron-Exon Ratio





Lab metrics

- Generated on-instrument, manually evaluated based on experience

Bioinformatic metrics

- From a population of libraries (minimum 50 runs)
 - determine 95th (warning) and 99th (fail) percentile





Metric


- A measured or calculated characteristic of a library

Threshold

- A value at which a library is to be assessed for quality

Not all metrics have thresholds

- Metrics that do not thresholds:
 - Read count; expected spike-in observed
- Metrics that have thresholds:
 - Reagent leftover, contamination rate, alignment rate



Hard Threshold


- Absolute point at which a library must be failed
- Indicates something has gone severely wrong
- Examples:
 - Very low alignment rate (<60%)
 - Very high contamination (>50%)

Outliers

- Metric beyond the 95th percentile of historical BCGSC data
- Contains usable data, but less than ideal
- Examples:
 - Low quality/low input material
 - Slightly lower genomic coverage
- BCGSC will manually review every library with 3 outlier metrics


**Why is QC
Important?**





How QC is useful to your processes

1. Confirm sample identity
 - Swaps or contamination events
2. Detect problems with laboratory processes
 - Uneven pooling, high ribosomal RNA content
3. To make improvements to protocols
 - How does a new protocol compare to the old version?
4. To compare results to previous experiments
 - Batch effects over time
 - Are additional lanes needed? How many?
5. To reduce costs
 - Avoid analyzing bad data and integrating results into existing data



How QC is useful to your science

- As a QC gate
 - Prevent bad data from being incorporated into an analysis
 - Sample swaps, low library diversity
- To identify outliers
 - Samples that have known issues that may affect analysis results
 - Explains observations in data when publishing results
- To perform trend analysis
 - Look at results over time
 - Provides a baseline by experiment type for comparison
 - Identify areas of optimization in lab & bioinformatic pipelines

Examples





Expected Indices	Observed Indices
TCCCGA	22%
ATCACG	26%
CTAGCT	24%
TGACCA	0%
No match	28%

Example 1:

- Conclusion – Incorrect 4th index specified
- Additional analysis – Examine no match bin
 - Infer missing index sequence based from most frequently observed index



Simple Examples: Index Splitting




Expected Indices	Library A Lane #1	Library A Lane #2	Library A Lane #3
TCCCGA	24%	0%	23%
ATCACG	20%	3%	18%
CTAGCT	23%	35%	28%
TGACCA	33%	10%	31%
No match	0%	52%	0%

Example 2:

- Conclusion – Lane #2 has been swapped with some other lane



Biological Exceptions – failing QC but the data is still usable



FFPE Samples

- Degraded DNA means PCR amplification was needed
- Higher duplicate read rate

Amplicon Libraries

- If amplicon sizes are small, high amounts of adapter are detected via read-through of fragment

Metagenomic Studies & Xenograft Libraries


- Alignment rate to a single target species may be low, but doesn't mean the data is bad

Low Input Libraries

- Frequently see higher background, lower fragment diversity



Safety Checks – when failing QC is a reason to stop



Low alignment rate

- BWA-aln works poorly on reads >125bp, use BWA-MEM
- Aligned to the wrong reference genome

Sample swaps

- Don't want to publish/analyze data for the wrong sample

Low bisulfite conversion rate of lambda phage


- Conversion reaction not done completely in lab

Genomic Contamination

- RNAseq library contains too much genomic DNA
 - Might affect observed expression levels

Conclusion





What is Quality Control

- 3 levels of QC

How QC is Carried out at the BCGSC


- Laboratory
- Automated Bioinformatic QC Pipeline
 - Role of manual review
- Some data that fails QC can sometimes be used

How is QC Useful

- Saves time in data analysis
- Aids in interpretation of data (publication)
- Identifies trends and areas for improvement



Acknowledgements



Genome Sciences Centre

Dr. Marco Marra

Dr. Steven Jones

Dr. Yussanne Ma

Bioinformatics Quality Control

- Eric Chuah
- Irene Li
- Gina Choe
- Dorothy Cheung
- Correy Lim
- Robert Lin

Laboratory Production Teams

- Dr. Andy Mungall
- Dr. Richard Moore
- Michael Mayo

GSC Production Teams

- Library Construction
- Sequencing Group
- LIMS
- BioApps Team
- Software Analysis
- Analysis Pipelines
- Data Analysis
- Systems Group
- Reanne Bowlby

More Information

BCGSC Website: <http://www.bcgsc.ca>

Email: ktse@bcgsc.ca