



computecanada
regional partner



UBC100

Advanced Research Computing

Education Outreach and Training Tutorials

Introduction to Short Read Mapping:
The foundation of next generation
sequencing analysis

April 3rd, 2019 (10:00AM-11:00PM PST)

Phillip A Richmond, Oriol Fornes

Copyright Information



Attribution-ShareAlike 3.0 Unported (CC BY-SA 3.0)

This is a human-readable summary of (and not a substitute for) the [license](#). [Disclaimer](#).

You are free to:

Share — copy and redistribute the material in any medium or format

Adapt — remix, transform, and build upon the material for any purpose, even commercially.

The licensor cannot revoke these freedoms as long as you follow the license terms.



The material is open source, and in this presentation no previous external work was utilized.



UBC100

Advanced Research Computing

Welcome!

- Welcome to the Introduction to Short Read Mapping
- I am co-teaching this seminar with Dr. Oriol Fornes, who studies gene regulation and frequently processes short-read data on the Cedar compute cluster.
- This is not meant to be a follow-along seminar, but the commands, datasets, and scripts will be available afterwards for your own exploration
- This presentation will be recorded and the slides will remain available

<http://bit.ly/2WD1ORc>



UBC100

Interactive Experience

We hope this is an interactive experience for all of you.

Questions/Problems can be posted to the Etherpad:

https://etherpad.openstack.org/p/EOT_APRIL2019

Dr. Oriol Fornes will be here to help answer questions while I'm presenting.



UBC100

Speaker Bio

Phillip Richmond

PhD Candidate, Wasserman Lab, BC Children's Hospital Research Institute

Bioinformatics Program, University of British Columbia

<https://phillip-a-richmond.github.io>

Research: Maximizing the Utility of Whole Genome Sequencing in the Diagnosis of Rare Genetic Disorders

Previous work in Genomics: Genomic Contributions to Ethanol Sensitivity in Mice, Polyploid Evolution in Yeast, Brewing Yeast Genomics, Cancer Cell Epigenetics, Addiction Predisposition

Also loves teaching genomics, and my puppy Sherlock Holmes



Session Outline


- Introduction to next generation sequencing data & diverse data types
- Transcriptional regulatory datasets and where to find them
- Mapping reads to the genome using BWA mem
- Peak calling and creating pileup files using MACS2
- Visualizing data in IGV

Session Outline


- Introduction to next generation sequencing data & diverse data types
- Transcriptional regulatory datasets and where to find them
- Mapping reads to the genome using BWA mem
- Peak calling and creating pileup files using MACS2
- Visualizing data in IGV

Next generation sequencing: Short-read sequencing

Fragments of DNA



Next generation sequencing: Short-read sequencing



Next generation sequencing: Short-read sequencing

Fragments of DNA




Adapter
Ligation

Sequencing Library




Sequencing
Reaction



1-Ligate to flowcell

Next generation sequencing: Short-read sequencing

Fragments of DNA





Adapter
Ligation


Sequencing Library



Sequencing
Reaction



1-Ligate to flowcell



2-Cluster amplify

Next generation sequencing: Short-read sequencing

Fragments of DNA




Adapter
Ligation


Sequencing Library




Sequencing
Reaction




1-Ligate to flowcell



2-Cluster amplify



3-Sequencing by Synthesis



Next generation sequencing: Short-read sequencing

Fragments of DNA




Adapter
Ligation


Sequencing Library




Sequencing
Reaction




1-Ligate to flowcell



2-Cluster amplify



3-Sequencing by Synthesis



Next generation sequencing: Short-read sequencing

Fragments of DNA





Adapter
Ligation


Sequencing Library




Sequencing
Reaction




1-Ligate to flowcell



2-Cluster amplify



3-Sequencing by Synthesis



Next generation sequencing: Short-read sequencing

Fragments of DNA





Adapter
Ligation

Sequencing Library



Sequencing
Reaction



1-Ligate to flowcell

2-Cluster amplify

3-Sequencing by Synthesis



```
@Read1  
TCTTGCGTACGTCTTCGATCGTA  
+  
!!@$@##@!%!@#$$!!LLBBDKSNK
```

Convert to
Fastq



UBC100

Advanced Research Computing

Diverse Input, Same Output Format

- Different inputs still result in the same output data format
- Examples:
 - DNA-seq, ChIP-seq, RNA-seq, GRO-seq, ATAC-seq
- For non-DNA assays (e.g. RNA-seq/GRO-seq), they undergo a conversion from RNA-->cDNA before sequencing

EXAMPLE

MEANING

```
@K00171:617:HMMTNBBXX:1:1101:28686:1648
1:N:0:GACTAGTA
TCTTGCGTACGTCTTCGATCGTA
+
!!@$@##@!%!@#$!LLBBDKSNK
```

```
@Readname:And:Flowcell:Info 1 or 2 for read pair:N:0:Barcode
Sequence
"Plus Sign"
ASCII-Quality Scores
```



Session Outline

- Introduction to next generation sequencing data & diverse data types
- Transcriptional regulatory datasets and where to find them
- Mapping reads to the genome using BWA mem
- Peak calling and creating pileup files using MACS2
- Visualizing data in IGV

ENCODE - Encyclopedia of DNA Elements

ENCODE is one of the many places to find open source data:

www.encodeproject.org



Based on an image by Darryl Leja (NHGRI), Ian Dunham (EBI), Michael Pazin (NHGRI)

You can download a diverse set of data across tissues/cell types

ENCODE Data Encyclopedia Materials & Methods Help

Experiment matrix

Enter search terms to filter the experiments included in the matrix.

Assay category

- DNA binding 8699
- Transcription 3132
- DNA accessibility 1100
- RNA binding 699
- DNA methylation 500

Assay

- ChIP-seq 8699
- DNase-seq 835
- polyA RNA-seq 770
- shRNA RNA-seq 523
- total RNA-seq 400

Experiment status

Selected filters: released

- released 14659
- archived 1044
- revoked 265

14659 results

Clear Filters

BIOSAMPLE

ASSAY

ChIP-seq DNase-seq polyA RNA-seq shRNA RNA-seq total RNA-seq eCLIP DNase array small RNA-seq WGBS microRNA-seq RNA microarray RAMPAGE RNA Bind-n-Seq genotyping array CAGE microRNA counts Replic-seq RBBS ...and 25 more

cell line	ChIP-seq	DNase-seq	polyA RNA-seq	shRNA RNA-seq	total RNA-seq	eCLIP	DNase array	small RNA-seq	WGBS	microRNA-seq	RNA microarray	RAMPAGE	RNA Bind-n-Seq	genotyping array	CAGE	microRNA counts	Replic-seq	RBBS
K562	669	10	19	268	11	245	3	7	1	2	10	1	2	9	1	50	6	1
HepG2	357	3	11	255	5	210	3	3	2	6	2	6	1	6	2	1	6	2
A549	374	14	27				2	9	1	5	2		2	3			2	1
GM12878	226	2	13	3			3	6	1	1	7	1	2	6	1		6	2
HEK293	257						2				1		2					2

+ See 254 more...

tissue

tissue	ChIP-seq	DNase-seq	polyA RNA-seq	shRNA RNA-seq	total RNA-seq	eCLIP	DNase array	small RNA-seq	WGBS	microRNA-seq	RNA microarray	RAMPAGE	RNA Bind-n-Seq	genotyping array	CAGE	microRNA counts	Replic-seq	RBBS
liver	162	14	20	3			1	1	9	7	7	2		7				1
stomach	106	21	15	5			3	4	10	4	6	5		4				1
heart	99	21	16	3			1	10	9	7	2			8				
lung	79	16	11	1			2	1	7	4	4	1		4				1
kidney	69	15	13				2	5	4	4				4				1

+ See 155 more...

whole organisms

whole organism	ChIP-seq	DNase-seq	polyA RNA-seq	shRNA RNA-seq	total RNA-seq	eCLIP	DNase array	small RNA-seq	WGBS	microRNA-seq	RNA microarray	RAMPAGE	RNA Bind-n-Seq	genotyping array	CAGE	microRNA counts	Replic-seq	RBBS
whole organism	1879	60	50											15				
carcass	8	4												4				

primary cell


primary cell	ChIP-seq	DNase-seq	polyA RNA-seq	shRNA RNA-seq	total RNA-seq	eCLIP	DNase array	small RNA-seq	WGBS	microRNA-seq	RNA microarray	RAMPAGE	RNA Bind-n-Seq	genotyping array	CAGE	microRNA counts	Replic-seq	RBBS
bone marrow-derived macrophage		14	78															
foreskin keratinocyte	37	2	5	3					3	13	13							

Chromatin Immunoprecipitation (ChIP-seq)

Chromatin Immunoprecipitation Sequencing (ChIP-seq) protocol:

Purpose: To find which sequences of DNA a specific protein interacts with, i.e Transcription Factor (TF).

1-Crosslink
DNA:Protein




Chromatin Immunoprecipitation (ChIP-seq)


Chromatin Immunoprecipitation Sequencing (ChIP-seq) protocol:

Purpose: To find which sequences of DNA a specific protein interacts with, i.e Transcription Factor (TF).

1-Crosslink
DNA:Protein



2-Shear




Chromatin Immunoprecipitation (ChIP-seq)


Chromatin Immunoprecipitation Sequencing (ChIP-seq) protocol:

Purpose: To find which sequences of DNA a specific protein interacts with, i.e Transcription Factor (TF).


1-Crosslink
DNA:Protein



2-Shear



3-Pull Down
protein using
anti-protein
antibody on a
column, wash
away other DNA




Chromatin Immunoprecipitation (ChIP-seq)


Chromatin Immunoprecipitation Sequencing (ChIP-seq) protocol:

Purpose: To find which sequences of DNA a specific protein interacts with, i.e Transcription Factor (TF).


1-Crosslink
DNA:Protein




2-Shear



3-Pull Down
protein using
anti-protein
antibody on a
column, wash
away other DNA



4-Reverse
Crosslink




Chromatin Immunoprecipitation (ChIP-seq)


Chromatin Immunoprecipitation Sequencing (ChIP-seq) protocol:

Purpose: To find which sequences of DNA a specific protein interacts with, i.e Transcription Factor (TF).


1-Crosslink
DNA:Protein




2-Shear




3-Pull Down
protein using
anti-protein
antibody on a
column, wash
away other DNA



4-Reverse
Crosslink



5-Ligate
sequencing
adapters




Chromatin Immunoprecipitation (ChIP-seq)


Chromatin Immunoprecipitation Sequencing (ChIP-seq) protocol:

Purpose: To find which sequences of DNA a specific protein interacts with, i.e Transcription Factor (TF).


1-Crosslink
DNA:Protein




2-Shear



3-Pull Down
protein using
anti-protein
antibody on a
column, wash
away other DNA




4-Reverse
Crosslink




5-Ligate
sequencing
adapters



6-Sequence
Library



ATAC-seq represents open chromatin



Nature Reviews | Genetics

Mapping data to a reference: ChIP-seq Peak Calling

- Individually, the short sequencing reads do not have much information
- Collectively, they can represent something useful
- Analyzing short-read data takes two common forms:
 - **Reference-based mapping**
 - Assembly

Example: ChIP-seq for a Transcription Factor




Raw data (not that useful)

Mapping data to a reference: ChIP-seq Peak Calling

- Individually, the short sequencing reads do not have much information
- Collectively, they can represent something useful
- Analyzing short-read data takes two common forms:
 - **Reference-based mapping**
 - Assembly


Example: ChIP-seq for a Transcription Factor



Mapping data to a reference: ChIP-seq Peak Calling

- Individually, the short sequencing reads do not have much information
- Collectively, they can represent something useful
- Analyzing short-read data takes two common forms:
 - Reference-based mapping
 - Assembly

Example: ChIP-seq for a Transcription Factor



Session Outline

- Introduction to next generation sequencing data & diverse data types
- Transcriptional regulatory datasets and where to find them
- Mapping reads to the genome using BWA mem
- Peak calling and creating pileup files using MACS2
- Visualizing data in IGV

For those of you who want to “follow along”

I'll be showing you files and pipeline scripts which are available to you to reuse/repurpose as you see fit.

NOTE: I do not expect you to follow along on the command line exploring the files. You can always re-watch this recording and hit pause. Maybe listen to it in slow motion. Whatever floats your boat.

Logging into cedar:

```
$ ssh username@cedar.computecanada.ca
```



UBC100

A place of learning

The main directory for today's workshop data/scripts is:

```
/scratch/richmonp/TRAINING/APRIL2019/
```

If you want, you can make a temporary directory here to play around with. If you do, name it something unique.

```
$ mkdir /scratch/richmonp/TRAINING/APRIL2019/SHERLOCK/
```

Change SHERLOCK to your own directory name if you want to rerun this script.

All you need is scripts

/scratch/richmonp/TRAINING/APRIL2019/SCRIPTS/ has 3 scripts inside it:

H3K27Ac_Workshop.sh

POLR2A_Workshop.sh

ATAC-Seq_Workshop.sh

I'm going to copy these so I can play with them:

```
$ cp /scratch/richmonp/TRAINING/APRIL2019/SCRIPTS/*sh SHERLOCK/
```

There are also scripts in this directory without the _Workshop, they are the ones I've already edited to work for my personal directory.

Breakdown of the script: Welcome to the mellow yellow

```
#!/bin/bash

#SBATCH --account=rrg-wyeth

## Mail Options
#SBATCH --mail-user=prichmond@cmmt.ubc.ca
#SBATCH --mail-type=ALL

## CPU Usage
#SBATCH --mem=50G
#SBATCH --cpus-per-task=16
#SBATCH --time=4-0:00
#SBATCH --nodes=1

## Output and Stderr
#SBATCH --output=%x-%j.out
#SBATCH --error=%x-%j.error
```

This header information contains info about the account to bill for these hours, I want it to mail me, how much RAM and CPUs I need over a single node, and where to send standard error and output

Breakdown of the script: Welcome to the mellow yellow

```
#!/bin/bash
```

```
#SBATCH --account=rrg-wyeth
```

```
## Mail Options
```

```
#SBATCH --mail-user=prichmond@cmmt.ubc.ca
```

```
#SBATCH --mail-type=ALL
```

```
## CPU Usage
```

```
#SBATCH --mem=50G
```

```
#SBATCH --cpus-per-task=16
```

```
#SBATCH --time=4-0:00
```

```
#SBATCH --nodes=1
```

```
## Output and Stderr
```

```
#SBATCH --output=%x-%j.out
```

```
#SBATCH --error=%x-%j.error
```

You will need to change these
to be relevant to your own use
case

This header information contains info about the account to bill for these hours, I want it to mail me, how much RAM and CPUs I need over a single node, and where to send standard error and output

Load my necessary tools

```
# Requirements
module load python/2.7.14
pip install numpy --user
pip install macs2 --user
module load bwa
module load samtools
MACS2=$HOME/.local/bin/macs2
```

I'm also going to load the necessary modules, and install a local version of MACS2 to my home directory.

Then, I set the MACS2 variable (blue guy) to be the command which calls the MACS2 tool. You'll see why later

You're going to need a reference genome next

```
# Genome
GENOME="/cvmfs/ref.mugqic/genomes/species/Homo_sapiens.GRCh38/genome/Homo_sapiens.GRCh38.fa"
INDEX="/cvmfs/ref.mugqic/genomes/species/Homo_sapiens.GRCh38/genome/bwa_index/Homo_sapiens.GRCh38.fa"
```

Next, I specify the genome I want to use to map my data against. I realize you won't all work in human, but if you work in a model organism I recommend checking out this repository for genomes:

`/cvmfs/ref.mugqic/genomes/`

Here, I'm using their **BWA index**, and their **Fasta file**

Reference Genome, Fasta file format

Reference genomes are packaged into fasta files.

Format:

```
>chromosome1_Name OtherChromInfo AccessionInfo Etc.
```

```
NNNNNNATTCTTGGATGGATAGCATGATCAGTAGACATGACATGACAGATGAGGGATATGATGACCA  
CCACCCAGATTCCCGGCCGGCCGGCCGGCCGGCCGGCCGGCCGGCCGGCCGGCCGGCCGGCCGGCTATATATATA  
CATAG ....
```

```
>chromosome2_Name OtherChromInfo AccessionInfo Etc.
```

```
NNNNNNNCCCCGGCCGGCCGGCCGGCCGGCCGGCCGGCCGGCCGGCCGGCCGGCCGGCCGGCTATATATATACAT  
AGATGATCAGTAGACATGACATGACAGATGAGGGATATGATGACCACCACCCAGATTGGAGTTGCCA  
GAT
```

We need to “index” this genome in order to map to it. There are many different genome indexing strategies. For `bwa`, we use the command `bwa index`, which creates an FM-Index of the genome.

```
$ bwa index <in.fasta>
```

This will generate these files:

```
genome.fa.amb, genome.fa.ann, genome.fa.bwt, genome.fa.pac, genome.fa.sa
```

Set some more variables

```
# Globals
ID="ATAC-seq"
WORKDIR=/scratch/richmonp/TRAINING/APRIL2019/DATA2/$ID
PLATFORM="Illumina"
SAMPLE="heart_left_ventricle"
THREADS=16
mkdir $WORKDIR
cd $WORKDIR
```

I'm setting an identifier, a working directory (change this if you want to use the script yourself), the sample name, and the threads I'm using. I also make a working directory and change into it.

I HIGHLY RECOMMEND using variables like this within your scripts. It will make it possible to easily change out a single variable or path, and the script can remain functional

Set some more variables

```
# Globals
ID="ATAC-seq"
WORKDIR=/scratch/richmonp/TRAINING/APRIL2019/DATA2/$ID
PLATFORM="Illumina"
SAMPLE="heart_left_ventricle"
THREADS=16
mkdir $WORKDIR
cd $WORKDIR
```

You will need to change this directory to be relevant to your own use case

I'm setting an identifier, a working directory (change this if you want to use the script yourself), the sample name, and the threads I'm using. I also make a working directory and change into it.

I HIGHLY RECOMMEND using variables like this within your scripts. It will make it possible to easily change out a single variable or path, and the script can remain functional

Set even more variables, and download some data

```
# Files
FASTQ_R1=$WORKDIR/$SAMPLE.R1.fastq.gz
FASTQ_R2=$WORKDIR/$SAMPLE.R2.fastq.gz
SAM_FILE=$WORKDIR/$SAMPLE.sam
BAM_FILE=$WORKDIR/$SAMPLE.bam
MACS2_DIR=$WORKDIR/MACS2
PEAKS_FILE=$MACS2_DIR/${SAMPLE}_peaks.narrowPeak

# Download ENCODE data
wget https://www.encodeproject.org/files/ENCFF766IGD/@download/ENCFF766IGD.fastq.gz
wget https://www.encodeproject.org/files/ENCFF075UOA/@download/ENCFF075UOA.fastq.gz
mv $WORKDIR/ENCFF766IGD.fastq.gz $FASTQ_R1
mv $WORKDIR/ENCFF075UOA.fastq.gz $FASTQ_R2
```

Here I set some file names, including for files that don't exist yet.

Then I download some data, and rename it according to the files I want them to be called.

If you want to explore lots of these datasets to download, use the www.encodeproject.org website.

Let the games begin: Mapping Reads to the Genome

```
# Map reads to genome
if [ ! -f $SAM_FILE ]; then
    # bwa mem [options] <idxbase> <in1.fq> [in2.fq]
    bwa mem -t $THREADS -M \
    -R "@RG\tID:$ID\tSM:$SAMPLE\tPL:$PLATFORM" \
    $INDEX $FASTQ_R1 $FASTQ_R2 > $SAM_FILE
fi
```

The little if/fi statements are to check if the output file exists, and if it does not exist, then perform the little command inside the block.

The BWA mem command is in the block, and at a minimum it needs an indexed genome, and an input fastq. I also add options -t for multithreading (using more cores), -R for a readgroup identifier (required for many tools), and -M for mapping split/secondary hits (not always needed). I also capture the standard out and place it into a SAM file.

The output SAM file

@SQ - Sequence (contig/chromosome) from reference file

@PG - Program information about mapping

@RG - Read group information (we won't have any here)

Tab delimited, each line is 1 read. Pairs will be next to each other in the file (e.g.

Line1: Read1

Line2: Read2

Col	Field	Type	Regex/Range	Brief description
1	QNAME	String	[!-?A-~]{1,254}	Query template NAME
2	FLAG	Int	[0,2 ¹⁶ -1]	bitwise FLAG
3	RNAME	String	* [!-()+-<>-~] [!-~]*	Reference sequence NAME
4	POS	Int	[0,2 ³¹ -1]	1-based leftmost mapping POSITION
5	MAPQ	Int	[0,2 ⁸ -1]	MAPping Quality
6	CIGAR	String	* ([0-9]+[MIDNSHPX=])+	CIGAR string
7	RNEXT	String	* = [!-()+-<>-~] [!-~]*	Ref. name of the mate/next read
8	PNEXT	Int	[0,2 ³¹ -1]	Position of the mate/next read
9	TLEN	Int	[-2 ³¹ +1,2 ³¹ -1]	observed Template LENGTH
10	SEQ	String	* [A-Za-z=.]+	segment SEQUENCE
11	QUAL	String	[!-~]+	ASCII of Phred-scaled base QUALity+33

Then we convert, sort, and index the bam file

```
# Convert SAM to sorted BAM
# create sorted BAM
if [ ! -f $BAM_FILE ]; then
    samtools view -Shu --threads `expr $THREADS - 1` $SAM_FILE | \
    samtools sort --threads `expr $THREADS - 1` -o $BAM_FILE -
    samtools index $BAM_FILE
fi
```

Here, I'm using the | to skip the step of saving the bam file, and then sorting it.

I link the two commands together to first convert the sam into bam using samtools view, and then sorting it using samtools sort.

I also add a multi-threading option, but samtools asks for “additional threads” so I take my thread# - 1.

The index command will create a .bai file next to the .bam file (file.bam.bai), which is needed for downstream tools

An easier version of samtools can be found here

```
$ module load samtools/1.3.1
```

We will use 3 samtools operations: view, sort, and index (in that order)

```
$ samtools view -b <in.sam> -o <out.bam>
```

```
$ samtools view -b Sample1.sam -o Sample1.bam
```

```
$ samtools sort <in.bam> -o <out.sorted.bam>
```

```
$ samtools sort Sample1.bam -o Sample1.sorted.bam
```

```
$ samtools index <in.sorted.bam>
```

```
$ samtools index Sample1.sorted.bam
```

Session Outline

- Introduction to next generation sequencing data & diverse data types
- Transcriptional regulatory datasets and where to find them
- Mapping reads to the genome using BWA mem
- Peak calling and creating pileup files using MACS2
- Visualizing data in IGV

The last component of the pipeline is to call peaks

```
# Call peaks
if [ ! -f $PEAKS_FILE ]; then
    $MACS2 callpeak -t $BAM_FILE \
        -n $SAMPLE \
        -p 0.00001 \
        --outdir $MACS2_DIR \
        -B --SPMR
fi
```

Here I'm calling peaks using MACS2.

I'm adding a sample name, I want it to output a bedgraph of normalized coverage for visualization. I'm using the ENCODE standard cutoff.

For ATAC-seq there is no "control", but for CHIP-seq pipelines there is sometimes a control sample, which you can provide as background for the peak caller.

Output files from MACS2

```
drwxrwxr-x 3 richmonp richmonp 4.0K Apr  2 22:31 ..
drwxrwxr-x 2 richmonp richmonp 4.0K Apr  2 22:46 .
-rw-rw-r-- 1 richmonp richmonp  87K Apr  2 22:51 NA_model.r
-rw-rw-r-- 1 richmonp richmonp 1.3G Apr  2 22:56 NA_treat_pileup.bdg
-rw-rw-r-- 1 richmonp richmonp 535M Apr  2 22:56 NA_control_lambda.bdg
-rw-rw-r-- 1 richmonp richmonp 1.2M Apr  2 22:56 NA_summits.bed
-rw-rw-r-- 1 richmonp richmonp 2.0M Apr  2 22:56 NA_peaks.xls
-rw-rw-r-- 1 richmonp richmonp 1.8M Apr  2 22:56 NA_peaks.narrowPeak
```

You'll get a set of files, and the ones which we will visualize are the:

*_treat_pileup.bdg,

*_summits.bed,

*_peaks.narrowPeak, which we will convert into a bed file (just rename it .bed)

Now, for visualization

I like to use OSX-Fuse / sshfs to connect my computer to Cedar. If you don't have it installed, google how to do so. There is also another 2-minute learn-along describing this process. If you ask me I'll dig it up.

```
Phillips-MacBook-Pro:~ philliprichmond$ sshfs richmonp@cedar.computecanada.ca:/home/richmonp/scratch/TRAINING/APRIL2019/DATA2/ ./Portal/  
richmonp@cedar.computecanada.ca's password:
```




UBC100

Then I'll open IGV


If you don't have IGV, I recommend downloading it here after this webinar:

<http://software.broadinstitute.org/software/igv/download>

Select the hg38 genome, if it isn't in your list..




Then go and get it



And then load your files via File> Load from file...

And you're going to want to select the OSXFUSE Volume 0 (sshfs)



Load in the .bam files, and the .narrowPeak files

The image shows a Mac desktop environment. In the background, the IGV (Integrative Genomics Viewer) application is open, displaying a track for 'Human (hg38)' with a search bar set to 'All'. In the foreground, a file browser window titled 'OSXFUSE Volume 0 (sshfs)' is open, showing a directory structure. The 'ATAC-seq' folder is expanded, listing several files. The 'NA_peaks.narrowPeak' file is highlighted in blue.

Name	Size
ATAC-seq	--
heart_left_ventricle.bam	4.43 GB
heart_left_ventricle.bam.bai	2.4 MB
heart_left_ventricle.fastq.gz_R1	2.07 GB
heart_left_ventricle.fastq.gz_R2	2.33 GB
heart_left_ventricle.sam	29.9 GB
MACS2	--
NA_control_lambda.bdg	560.3 MB
NA_model.r	88 KB
NA_peaks.narrowPeak	1.8 MB

And explore away!

This is a good region for the heart transcriptional regulation:

Chr10:21160000-23400000

Around the gene NEBL.

I'll now explore this data interactively and open to questions / comments.

To visualize the bedgraph files effectively...

You'll need to convert them to bigwig

This can be done, and if you're interested in learning how let me know.

In fact, Oriol has been working on it this morning so we should have the command ready for you soon!

:)

Open question and answer period



UBC100

Acknowledgements

- Phil Richmond (Teacher)
 - PhD Student Wasserman Lab, enjoys teaching
- Oriol Fornes (Co-teacher)
 - Post-doc, Deputy Group Leader, Wasserman Lab

FLASH DEBUGGING

```
$ samtools sort Sample1.bam -o Sample1.sorted.bam
```

Crazy characters printing to the screen

```
$ samtools view -bS Sample1.sam Sample1.bam
```

Crazy characters printing to the screen

```
$ samtools index Sample1.bam
```

[E::hts_idx_push] unsorted positions

samtools index: "Sample1.bam" is corrupted or unsorted

```
$ bwa mem -t ../GENOME/genome.fa Sample_R1.fastq
```

```
Sample_R2.fastq
```

[E::bwa_idx_load_from_disk] fail to locate the index files

Fix: This sort command doesn't use a -o

Unless you specify -T and -O as well.

```
$ samtools sort Sample1.bam Sample1.sorted
```

Fix: This commands needs a -o for the output

```
$ samtools view -bS Sample1.sam -o Sample1.bam
```

Fix: Order matters. Sort before you index

```
$ samtools index Sample1.sorted.bam
```

Fix: the -t option requires an integer. Otherwise, all the other positional arguments are out of place.

```
$ bwa mem -t 4 ../GENOME/genome.fa Sample_R1.fastq  
Sample_R2.fastq
```



ERROR: Loading SAM/BAM index files are not supported: /Users/philliprichmond/Desktop/NA20845.chr19.subregion.sorted.bam.bai
Load the SAM or BAM file directly.



Fix: Make sure you load the .bam file,
The .bai file just needs to be in the same directory
As the .bam file



UBC100

Advanced Research Computing